
# DISPERSA

## HPMC VISCOELASTIC

### Rheological Profile of DISPERSA

Measured at 25°C

## DESCRIPTION

DISPERSA HPMC Viscoelastic is a sterile, high molecular weight, hydroxypropyl methylcellulose solution for use as a surgical aid in anterior segment surgical procedures including cataract extraction and intraocular lens implantation. DISPERSA HPMC Viscoelastic protects the corneal endothelial cells during ophthalmic surgical procedures and maintains a deep anterior chamber. DISPERSA HPMC Viscoelastic is intended to be removed from the eye and can be easily irrigated/aspirated from the anterior chamber following surgery. In addition, DISPERSA HPMC Viscoelastic can be used to coat surgical instruments for added protection of ocular tissues and endothelium cells. DISPERSA is optically clear and can be stored at room temperature.

DISPERSA HPMC Viscoelastic is available in viscosities of 8,000 cps and 20,000 cps (shear rate 0.1/second). It is supplied in disposable 3cc syringes along with a with a 23-gauge (thin wall) stainless steel delivery cannula. The syringes are medical grade

polypropylene. The syringe has medical grade elastomeric stoppers.

Each 1 mL of DISPERSA HPMC Viscoelastic contains approximately 16 to 20 mg/mL hydroxypropyl methylcellulose (HPMC) premium USP/EP dissolved in a physiologically compatible buffering solution (BSS) containing 0.48 mg Calcium Chloride Dihydrate USP, 0.30 mg Magnesium Chloride Hexahydrate USP, 0.75 mg Potassium Chloride USP, 3.90 mg Sodium Acetate Trihydrate, 6.40 mg Sodium Chloride USP, and 1.70 mg Sodium Citrate Dihydrate USP in Water for Injection USP. The final product pH is 6.8 to 7.6 and the Osmolality is 200 to 400 milliOsmoles/kilogram.

## INTENDED USE

DISPERSA HPMC viscoelastic is intended for use as a surgical aid in anterior segment surgical procedures including cataract extraction and intraocular lens implantation. DISPERSA HPMC Viscoelastic protects the corneal endothelial cells during ophthalmic surgical procedures and maintains a deep anterior chamber. DISPERSA HPMC Viscoelastic is intended to be removed from the eye and can be easily irrigated/aspirated from the anterior chamber following surgery. In addition, DISPERSA can be used to coat surgical instruments and intraocular lenses prior to insertion into the eye for added protection of ocular tissues and endothelium cells.

## CONTRAINDICATIONS

Presently, there are no known contraindications to the intraocular use of hydroxypropyl methylcellulose viscoelastic when used as recommended.

## PRECAUTIONS

Precautions associated with the use of DISPERSA relate to the surgical procedure being performed. The following precautions are recommended:

1. Do not over fill the anterior chamber of the eye with DISPERSA. Over filling may cause increased intraocular pressure, glaucoma, or other associated ocular damage.

2. Remove as much as possible of the DISPERSA at the close of surgery by irrigation/aspiration. Use care to avoid injury to the corneal endothelial cells as a result of the removal procedure.

3. Monitor postoperative intraocular pressure. If significant rises occur, treat with appropriate therapy. Postoperative intraocular pressure should be evaluated as a result of pre-existing glaucoma, impeded outflow, and as a result of related operative procedures which may include enzymatic zonulyses, absence of an iridectomy, trauma to filtration structures, and by blood and lenticular remnants in the anterior chamber.

4. The concurrent presence of medication(s) in the anterior chamber or associated ocular structures may interact with DISPERSA to cause clouding.

5. DISPERSA is preservative free and should be discarded after a single use. Avoid using any excess viscoelastic. Do not resterilize.

6. DISPERSA should be used at ambient temperature. Protect from direct light and freezing and excessive heat greater than 45°C.

7. DO NOT REUSE. This is a single use, disposable, device that cannot be properly cleaned and re-sterilized.

## ADVERSE REACTIONS

Hydroxypropyl methylcellulose viscoelastic is tolerated well after injection into human eyes. A transient rise of intraocular pressure postoperatively has been reported. Isolated incidences of postoperative inflammatory reactions (iritis, hypopyon) as well as corneal edema and corneal decompensation have also been reported. Their

relationship to hydroxypropyl methylcellulose viscoelastics has not been established.

## CLINICAL APPLICATIONS

Remove the container of DISPERSA from the pouch. DISPERSA is guaranteed to be sterile unless the pouch is damaged or opened. A delivery cannula is also packaged in the pouch. Attach the cannula to the syringe. Bleed the air from the cannula by ejecting a small amount of viscoelastic from the cannula.

Slowly introduce a sufficient amount of DISPERSA into the anterior chamber to fill the eye. The injection of DISPERSA can be performed either before or after removal of the crystalline lens. Injection of DISPERSA prior to lens removal will have the additional advantage of protecting the corneal endothelium from possible damage during the removal process. Additional volumes of DISPERSA can be injected during the procedure to replace any viscoelastic lost during surgery.

DISPERSA may also be used to coat the surgical instruments and the intraocular lens prior to insertion into the eye. It may also be used to lubricate instruments such as folded lens inserters.

## SAFETY & EFFECTIVENESS

Extensive pre-clinical and clinical testing has been performed on hydroxypropyl methylcellulose viscoelastics. This material has been shown to be safe and effective for use as an intraocular viscoelastic. DISPERSA HPMC viscoelastic has been tested and shown to be sterile, non-pyrogenic, non-cytotoxic, non-irritating in the rabbit anterior chamber irritation test, and free of potentially antigenic proteinaceous material.

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